GE Healthcare

Serial Interface Data Services Service Manual

Serial Interface Data Services
English
© 2005 General Electric Company
All Rights Reserved.

NOTE: Due to continuing product innovation, specifications in this manual are subject to change without notice.

Listed below are GE Medical Systems *Information Technologies* trademarks. All other trademarks contained herein are the property of their respective owners.

DASH, SOLAR, and UNITY NETWORK are trademarks of GE Medical Systems *Information Technologies* registered in the United States Patent and Trademark Office.

UNITY is a trademark of GE Medical Systems *Information Technologies*.

Contents

1

Introduction 1-1

Manual Information	1-2
Intended Use	1-2
Ordering Manuals	1-2
Revision History	1-2

2

Hardware Connections 2-1

Overview	2-2
----------------	-----

Connection Specifications	2-2
---------------------------------	-----

Data Retrieval	2-3
Serial Port Limitations	2-3
Data Request Packet	2-3
Request Packet Definition	2-4
Functional Specifications	2-5
Function and Subfunction Codes	2-6
Response Packet Definition	2-7

3

Parameters 3-1

Client/Server Communications Model	3-2
Defining the Client/Server System	3-2
Queries/Responses	3-2

Data Packet Architecture	3-3
Parameter Update Packet	3-3
Data Representation	3-4
Reserved Data Values	3-4
Bedside Message Structure	3-5
Bedside Float Structure	3-6
Parameter Float Structure	3-7
Parameter Update Structure	3-8
Extended Parameter Update Structure	3-8
Setup and Limits Structure	3-9
Message Structure	3-9
More Setup Structure	3-10
Miscellaneous Data Structures	3-10

ECG Parameter	3-11
Parameter Update	3-11
Extended Parameter Update	3-11
Setup and Limits	3-11
Messages	3-12
More Setup	3-13
Miscellaneous	3-13
ST Parameter	3-14
Parameter Update	3-14
Extended Parameter Update	3-14
Setup and Limits	3-14
Messages	3-14
More Setup	3-14
Miscellaneous	3-14
12-Lead ST Parameter	3-15
Parameter Update	3-15
Extended Parameter Update	3-15
Setup and Limits	3-15
Messages	3-16
More Setup	3-16
Miscellaneous	3-16
Respiration Parameter	3-17
Parameter Update	3-17
Extended Parameter Update	3-17
Setup and Limits	3-17
Messages	3-18
More Setup	3-18
Miscellaneous	3-18
Blood Pressure Parameter	3-19
Parameter Update	3-19
Extended Parameter Update	3-19
Setup and Limits	3-20
Messages	3-20
More Setup	3-21
Miscellaneous	3-21
SpO₂ Parameter	3-22
Parameter Update	3-22
Extended Parameter Update	3-22
Setup and Limits	3-22
Messages	3-23
More Setup	3-23
Miscellaneous	3-23
Temperature Parameter	3-24
Parameter Update	3-24
Extended Parameter Update	3-24

Setup and Limits	3-24
Messages	3-25
More Setup	3-25
Miscellaneous	3-25
Cardiac Output Parameter	3-26
Parameter Update	3-26
Extended Parameter Update	3-26
Setup and Limits	3-26
Messages	3-27
More Setup	3-27
Miscellaneous	3-27
Non-invasive Blood Pressure Parameter	3-28
Parameter Update	3-28
Extended Parameter Update	3-28
Setup and Limits	3-28
Messages	3-29
More Setup	3-29
Miscellaneous	3-29
CO₂ Parameter	3-30
Parameter Update	3-30
Extended Parameter Update	3-31
Setup and Limits	3-31
Messages	3-32
More Setup	3-33
Miscellaneous	3-33
SvO₂ Parameter	3-34
Parameter Update	3-34
Extended Parameter Update	3-34
Setup and Limits	3-34
Messages	3-35
More Setup	3-35
Miscellaneous	3-35
Ventilator Parameters	3-36
VENT_PAR Parameter	3-36
VENT_PAR1 Parameter	3-38
VENT_PAR2 Parameter	3-39
VENT_PAR3 Parameter	3-40
VENT_PAR4 Parameter	3-41
VENT_PAR5 Parameter	3-42
Gas Parameters	3-43
Parameter Update	3-44
Extended Parameter Update	3-45
Setup and Limits	3-45
Messages	3-46
More Setup	3-47

Miscellaneous	3-47
Arterial Blood Gas Parameter	3-48
Transcutaneous CO₂ / O₂ (Interfaced) Parameter	3-51
Parameter Update	3-51
Extended Parameter Update	3-51
Setup and Limits	3-51
Messages	3-52
More Setup	3-52
Miscellaneous	3-52
Continuous Cardiac Output (Interfaced) Parameter	3-53
Parameter Update	3-53
Extended Parameter Update	3-53
Setup and Limits	3-53
Messages	3-54
More Setup	3-54
Miscellaneous	3-54
IV Pump (Interfaced) Parameter	3-55
Parameter Update	3-55
Extended Parameter Update	3-55
Setup and Limits	3-55
Messages	3-55
More Setup	3-56
Miscellaneous	3-56
Urometer (Interfaced) Parameter	3-57
Parameter Update	3-57
Setup and Limits	3-57
Messages	3-57
More Setup	3-57
Miscellaneous	3-57
Pulse Oximeter (Interfaced) Parameter	3-58
Parameter Update	3-58
Extended Parameter Update	3-58
Setup and Limits	3-58
Messages	3-59
More Setup	3-59
Miscellaneous	3-59
ECG (Interfaced) Parameter	3-60
Parameter Update	3-60
Extended Parameter Update	3-60
Setup and Limits	3-60
Messages	3-60
More Setup	3-60
Miscellaneous	3-60
Blood Pressure (Interfaced) Parameter	3-61

Parameter Update	3-61
Extended Parameter Update	3-61
Setup and Limits	3-61
Messages	3-62
More Setup	3-62
Miscellaneous	3-62
Temperature (Interfaced) Parameter	3-63
Parameter Update	3-63
Extended Parameter Update	3-63
Setup and Limits	3-63
Messages	3-63
More Setup	3-63
Miscellaneous	3-63
NBP (Interfaced) Parameter	3-64
Parameter Update	3-64
Extended Parameter Update	3-64
Setup and Limits	3-64
Messages	3-64
More Setup	3-64
Miscellaneous	3-64
Respiration (Interfaced) Parameter	3-65
Parameter Update	3-65
Extended Parameter Update	3-65
Setup and Limits	3-65
Messages	3-65
More Setup	3-65
Miscellaneous	3-65
Blood Temperature/Cardiac Output (Interfaced) Parameter	3-66
Parameter Update	3-66
Extended Parameter Update	3-66
Setup and Limits	3-66
Messages	3-66
More Setup	3-66
Miscellaneous	3-66
Respiratory Mechanics Parameters	3-67
RM_PAR Parameter	3-67
RM_PAR1 Parameter	3-69
SvO₂ (Interfaced) Parameter	3-70
Parameter Update	3-70
Extended Parameter Update	3-70
Setup and Limits	3-70
Messages	3-70
More Setup	3-70
Miscellaneous	3-70

ICG Parameter	3-71
icg_par Parameter	3-71
icg1_par Parameter	3-73
icg2_par Parameter	3-74
BIS Module Parameter	3-75
Parameter Update	3-75
Extended Parameter Update	3-75
Setup and Limits	3-75
Messages	3-76
More Setup	3-76
Miscellaneous	3-76
EEG Parameter	3-77
eeg_par Parameter	3-77
eeg1_par Parameter	3-79
eeg2_par Parameter	3-80
eeg3_par Parameter	3-81
eeg4_par Parameter	3-82
eeg5_par Parameter	3-83

1 Introduction

Manual Information

Intended Use

This manual provides the serial interface data services information that an institution's information technology personnel can use to acquire parameter data from GE Medical Systems *Information Technologies* bedside monitors, including:

- Dash 2000/3000/4000/5000 Patient Monitor
- Solar 8000M Patient Monitor
- Solar 8000i Patient Monitor
- Solar 9500 Information Monitor
- Unity Network Interface Device

NOTE

Not all monitors support all of the parameters described in this manual. Refer to the monitor's operator's manual for more information on supported parameters.

Ordering Manuals

A paper copy of any manual will be provided upon request. Contact your local GE Medical Systems *Information Technologies* representative and request the part number on the first page of the manual.

Revision History

Each page of this document has the document part number and revision letter at the bottom of the page. The revision letter changes whenever the document is updated.

Revision	Comments
A	Initial release.

2 Hardware Connections

Overview

This chapter describes what is needed to physically connect a workstation or personal computer to any GE Medical Systems *Information Technologies* bedside monitor and how to retrieve data via the serial port.

Connection Specifications

There are three different types of serial port connectors. These communication ports all use 9600 baud rate, 8 data bits, one stop bit, and no parity.

		Solar 8000M/i	Solar 9500	Unity ID	Dash 2000/3000/4000/5000 Dash Port Docking Station
Standard		EIA RS-232			
Connector	Name	RS-232 1	10101 2 10101 1	RS 232	AUX, AUX 1, AUX 2
	Type	DB-9M			RJ-45
Required Pins	Transmit Pin	2	3	2	6
	Receive Pin	3	2	3	3
	RETURN	5			4
Isolation Provided		None provided. Electrical isolation is the responsibility of the 3rd party interfacing to the monitor.		Basic Insulation @ 250V	
Interconnect Cabling		Standard NULL Modem RS-232 cable with shield. 100 foot maximum length.			PC Interface DIDCA (PN 420915-013) with standard category 5 cable. 50 foot maximum length
Notes		Serial Data Services is only available on the RS-232 1 port. Refer to the Solar 8000M/i Service Manual for more information.	Serial Data Services can be configured for RS-232 port 1 or 2 but not both simultaneously. Refer to the "Processing Unit/ Polled Parameter Interconnection" section of the Solar 9500 Service Manual for more information.	Refer to the Unity Network ID Service Manual for more information.	Refer to the Dash Service Manual for detailed Auxiliary Communication Connector information. Refer to the Dash Port Docking Station or Dash Port 2 Docking Station service manual for additional information.


Data Retrieval

Serial Port Limitations

Data available from the serial port is limited to data from the monitor. Data from other monitors cannot be accessed even if the monitor is connected to a GE Medical Systems *Information Technologies* Unity Network. In addition, no waveform data is available through the serial port. The serial port does not have the necessary bandwidth to transmit the waveform data.

Data Request Packet

To acquire data from the serial port, a request packet must be sent to the monitor. The request packet is based upon the SBEDSIDE_MSG_DEF structure, which is defined in Chapter 3, “**Parameters**”. The following example shows a sample request packet:

Request Packet Definition

dst_addr	The first byte of the destination address is always set to 64 when using the serial port connection.
src_addr	The source address is not necessary. Set it to 0.
fun_code	The function code specifies what action the server is to perform. To simply read data from the server a function code of 202 would be used.
sub_code	The subfunction code further defines the request being sent to the server. In the above example the subfunction code 35 is sent to request polled parameters.
version	Determines the message structure to be used. For serial port communications this value is always set to 1.
seq_num	Not used, set to 0.
req_res	Not used, set to 0.
proc_id	Not used, set to 0.
oln	Not used, set to 0.
return_status	Not used, set to 0.
data_count	When a request is sent to the server the data count is set to 0.
CRC	Used to verify each received packet, to ensure data integrity.

Functional Specifications

Item	Description
Data format	Asynchronous serial
Transmission modes	Full duplex, half duplex
Transmission speed	9600 bits per second for Data Services connection
Packet structure	Data bits 8; Parity None; Stop bits 1; Speed 9600 bps
Error detection	CRC16
Polling frequency	Not more than once every 2 seconds
OSI model layers:	
Application	Parameter data structure
Presentation	BEDMSG structure
Transport	n/a
Network	n/a
Data link	RS-232 UART
Physical	Serial interface or standard category 5 cable

Function and Subfunction Codes

Function Codes

time broadcasts from time master

201	FC_WRITE	write request
202	FC_READ	read request
203	FC_ABORT	abort requests

Subfunction Codes

time broadcast update or write request

2	SC_ADMIT
3	SC_DISCHARGE
5	SC_PATIENT_ID
6	SC_PATIENT_NAME
10	SC_SOFTWARE_REVISION
35	SC_POLLED_PARAMETER_REQUEST

Response Packet Definition

The monitor responds to request packets quickly, returning a response packet. A response packet can contain the data that was requested or an acknowledgment of the request.

**Sample Response Packet for
Serial Interface Communication**

destination dst_addr	0	64	0	0	0	
	4	0	0	0	0	source src_addr
function fun_code	8	0	0	0	0	subfunction sub_code
version	12	0	201	0	20	sequence seq_num
req response req_res	16	0	1	0	0	process id proc_id
	20	0	0	0	0	
	24	0	0	0	0	
	28	0	0	0	0	
	32	0	0	0	0	origin locate oin
	36	0	0	0	0	
return_status	40	0	0	0	0	
	44	0	0	0	0	data count data_count
data field data	48	0	0	0	0	

NOTE

The first byte of the response packet coming from a Solar 9500 with version 4 software is 0x00, not 0x40.

The following elements are in a response packet:

- The first 60 bytes of the response packet follow the SBEDSIDE_MSG_DEF structure.
- The function code would be 201 which identifies a write command, instructing the requesting computer/workstation that the monitor is writing data to it.
- The subfunction code would be 20, specifying that a data stream containing parameter data is being sent.
- The version number is set to 1 to verify that the correct data structures are being used by the receiving computer/workstation. The other fields in the packet are not relevant and may be ignored until the data count field is received.
- The data count field specifies the quantity of data that follows in the data structures. Note that the value in data count refers only to the number of bytes in the data portion of the packet. It does not include the first 60 bytes of the packet, nor does it include the two-byte CRC that is added to the end of the packet.
- Following the data count field is the actual data array (data structures). The data array contains the requested patient data and alarms status information.
- The two-byte CRC value follows the data array.

3 Parameters

Client/Server Communications Model

This chapter describes the basic organization of data structures in the transmitted data packets.

Defining the Client/Server System

Each monitor acts as the server entity in a client/server environment. Your personal computer/workstation functions as the client, running the software application you program. In the following description the monitor is referred to as the server and your personal computer/workstation as the client.

Queries/Responses


Each time data requested by the client application, the application must issue an appropriate query to the server. This query is sent in the form of a request packet. In turn the server responds with data or at least an acknowledgment in the form of a response packet. Both the query and response are transmitted in the form of data packets. This section describes how these packets are constructed. Your client application must be written to construct similar query packets and parse received packets to recover the patient data.

Data Packet Architecture

This section provides a general description of the data communication packets used to construct and parse the serial interface data structures. Several structures are used to define the data communication and parameter packet architecture. Both the request and response packets are similarly constructed. In general, the response data packets will contain more data.

The packet architecture contains three types of nested parameter data structures. The Bedside Message Structure identifies the source and destination for the Bedside Float Structure, a function to be performed, and other specifications relating to the function and to the data that may be contained in the packet. The Bedside Float Structure holds the Parameter Float Structure with the Parameter Update packet data inside.

Parameter Update Packet

Each parameter (SPAR_FLOAT structure) contains data for a single parameter. The SPAR_FLOAT data structure is comprised of 5 data structures defined below. Use these data structure to access specific parameter and sub-parameter data. The illustration shows the next layer of nested structures that reside within the Bedside Message Structure. All these nested structures and subsequent structures within them are collectively referred to as the parameter update packet.

The Bedside Float Structure provides device level information, such as alarm state, alarm level, patient admission, and graph status. It also specifies the number of parameters that are included in the subsequent parameter data array

A Parameter Float Structure is present for each parameter in the Bedside Float Structure. A typical Bedside Float Structure will have many Parameter Float Structures contained within it. The Parameter Float Structure specifies the data structure for each parameter.

Each parameter float structure contains five substructures. Each substructure defines how to interpret the data for each parameter. Each parameter has different combinations of flags and data values. Definitions for the structures described above are provided later in this chapter.

For asynchronous communications, all response packets are followed by a 2-byte checksum. All received data should be checked by generating a CRC value and comparing it to the transmitted value.

Data Representation

The following data types are used throughout this document.

Description	Data Values
UTINY	unsigned 8-bit byte
CHAR	signed 8-bit byte
COUNT, SHORT	signed 16-bit word
UCOUNT	unsigned 16-bit word

Reserved Data Values

The following data values are used to represent unique parameter data conditions and should not be displayed or trended as actual physiologic data values.

```
#define INVALID      -32768      /* Invalid Data */
#define MISSING      -32767      /* Missing Data */
#define PAR_DRAW     -32766      /* Parameter Draw */
#define PAR_FLUSH    -32765      /* Parameter Flush */
#define PAR_ZERO     -32764      /* Parameter Zero */
#define PAR_CAL      -32763      /* Parameter Calibration */
#define NO_BP_PULSE  -32762      /* NBP no pulse */
#define SENSOR_FAIL  -32761      /* Parameter sensor fail */
#define INVALID_BYTE -128        /* Invalid Data (1 byte value) */
#define MISSING_BYTE -127        /* Missing Data (1 byte value) */
```

Bedside Message Structure

Bedside Message is parsed according to the SBEDSIDE_MSG_DEF structure to determine the destination, source, function to be performed, and the amount of data (assuming a response packet). This data is contained in a Bedside Float Structure.

```
typedef struct sbedside_msg_def
{
#define ASYNC_PDMS_FC          0x40          /* dst_addr[0] = ASYNC_PDMS_FC */
    UTINY dst_addr[6];          /* destination address */
    UTINY src_addr[6];          /* source address */
    COUNT fun_code;             /* function code */
#define FC_WRITE                201          /* write requests pat name etc */
#define FC_READ                 202          /* read software revisions etc */

    COUNT sub_code;             /* subfunction code */
    COUNT version;              /* version of bed_msg */
#define BEDMSG_CS_VER_5        0
#define BEDMSG_CS_VER_6        6

    COUNT seq_num;              /* response sequence number */
    COUNT req_res;              /* request response flag */
    COUNT proc_id;              /* requestors process id */
    UTINY oln[32];              /* origin location name */
    COUNT return_status;        /* return status */
    COUNT data_count;           /* following message data count */
    COUNT data[1];              /* parameter data(beginning SBEDSIDE_FLOAT) */
} SBS_MSG_DEF, *pSBS_MSG_DEF;
```

Bedside Float Structure

The Bedside Float Structure is parsed according to the SBEDSIDE_FLOAT structure to obtain additional device status data and to determine how many parameters are included in the subsequent data array.

```
typedef struct sbedside_float          /* length without par_float = 6 bytes */
{
    UTINY alarm_state;                /* active, silence, pause, off */
    #define ALARM_ACTIVE                0 /* alarms processed if patient admitted */
    #define ALARM_SILENCE                1 /* 1 minute graph, audio alarm hold */
    #define ALARM_PAUSE                2 /* 5 minute graph, audio alarm hold */
    #define ALARM_OFF                3 /* no alarms processed */
    #define ALARM_VOLUME_OFF            4 /* Make a sound every 3 minutes */
    #define ALARM_PAUSE_DISP_OFF        5 /* alarms paused, display off */

    UTINY alarm_level;                /* highest parameter alarm level */
    #define ALARM_LEVEL_STATUS_ONLY      0 /* no processing */
    #define ALARM_LEVEL_SYSTEM_MESSAGE  1 /* normal display, no audio */
    #define ALARM_LEVEL_SYSTEM_ADVISORY 2 /* fog_horn 1 tone audio alarm */
    #define ALARM_LEVEL_SYSTEM_WARNING  3 /* fog-horn continuous audio alarm */
    #define ALARM_LEVEL_MESSAGE          4 /* no audio */
    #define ALARM_LEVEL_ADVISORY         5 /* 1 beep non-latch audio alarm */
    #define ALARM_LEVEL_WARNING          6 /* 2 beep non-latch audio alarm */
    #define ALARM_LEVEL_CRISES           7 /* 3 beep latching audio alarm */

    UTINY audio_alarm_level;          /* alarms.h, current audio */

    UTINY patient_admission;          /* admitted or discharged */
    #define DISCHARGED                0
    #define ADMITTED                  1

    UTINY number_of_parameters;        /* length of par_float array */
    UTINY graph_status_msg;            /* not available for serial interface */
    SPAR_FLOAT par_float_list[1];      /* array of individual par float structures */
} SBEDSIDE_FLOAT, *pSBEDSIDE_FLOAT;
```

The amount of SPAR_FLOAT parameter data structures that must be parsed is determined from the *number_of_parameters* value in the SBEDSIDE_FLOAT structure.

NOTE

The maximum *number_of_parameters* is limited to 16.

Parameter Float Structure

Each parameter (SPAR_FLOAT structure) contains data for a single parameter. The SPAR_FLOAT data structure is comprised of 5 data structures defined below. Use these data structure to access specific parameter and sub-parameter data.

- Parameter update (PAR_UPD) structure,
- Extended parameter update (EXTENDED_PAR_UPD) structure,
- Setup and limits (LIMIT_VALUES) structure,
- More setup (MORE_SETUP) structure, and
- Messages (PAR_MSG) structure.

Each data structure contains two unsigned 8-bit values named PAR_FUNC_CODE and PARCODE.

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Parameter Update Structure

This structure is intended to provide the current parameter values and some basic status information.

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The PAR_STATUS field is typically a bit field (16 bits) which provides flags indicating limit violations, alarm status, and status of the data acquiring device. Definitions for these flags are provided.

The PAR_VALUE array holds the current parameter values. The number of PAR_VAL entries used depends upon the parameter. The PAR_FUNC_CODE for parameter update is 1.

Extended Parameter Update Structure

This structure is for those parameters which require additional locations for information.

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Setup and Limits Structure

This structure communicates setup information and limit values.

```
struct LIMIT_VALUES
{
    COUNT    lo_limit;
    COUNT    hi_limit;
};

struct SETUP_N_LIM
{
    UTINY                                par_func_code;
    UTINY                                parcode;
    UCOUNT                                flag[2];
    struct LIMIT_VALUES                  limit_values[3];
    COUNT                                extra_limit;
};
```

The flag array is used differently by each parameter to convey various setup data. The LIMIT_VALUES array typically provides limit values for corresponding parameter values. The number of LIMIT_VALUES entries used depends upon the parameter. The EXTRA_LIMIT field provides space for one extra limit value. The PAR_FUNC_CODE for parameter update is 3.

Message Structure

This structure communicates information specifying a specific message that should be displayed with reference to the given parameter. The messages are defined and identified with a message index/code. This message index/code is sent in the PAR_MSG array.

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More Setup Structure

This structure is for those parameters which require additional setup information beyond what is included in the Setup and Limits structure. The VAL array provides 8 additional bytes of data. This data is used differently by each parameter. The PAR_FUNC_CODE for more setup is 2.

```
struct MORE_SETUP
{
    UTINY    par_func_code;
    UTINY    parcode;
    COUNT    val[4];
};
```

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Messages

```

par_func_code (char)PAR_MSG_FC (21)
parcode (char) hr_par (58)
    messages (struct PAR_MSG [3])
        messages[0].msg_index
        messages[1].msg_index
            CODENONE      0      /* no message */
            CODENORM      1      /* Normal rhythm */
            CODELCFAV     2      /* Low count (<20) of favorite beat this min */
            CODEEPNC      3      /* Pace non-capture */
            CODEEPNS      4      /* Pace non-sense */
            CODEAFIB      5      /* Atrial fib -> irregular */
            CODEAFLT      6      /* Atrial flutter */
            CODEBRAD      7      /* Sinus brady */
            CODETACH      8      /* Sinus tach */
            CODEPSVC      9      /* Isolated PSVC */
            CODEPSVP     10      /* PSVC pair */
            CODESVB      11      /* SV brady */
            CODESVT      12      /* SV tach */
            CODESTDV     13      /* ST deviation */
            CODEPVC      14      /* Isolated PVC */
            CODETGMY     15      /* Trigeminy */
            CODEPAUS     16      /* Missing beat */
            CODEACCV     17      /* Accel. ventricular */
            CODEBGMY     18      /* Bigeminy */
            CODECPLT     19      /* Couplet */
            CODEVBRD     20      /* V brady */
            CODERONT     21      /* R on T */
            CODEVT35     22      /* Short run of V tach */
            CODEVTAC     23      /* V tach */
            CODEVFIB     24      /* V fib */
            CODEASYS     25      /* Asystole */
            CODEAFIB_EKPROV10 28 /* Atrial Fibrillation - using AFIB algorithm in
                                   EkPro v10.1 */

    /* Following are Telemetry Specific */
        CODENOTELEM     26      /* TTX no telem */
        CODENOTELEM_MSG_ONLY 27 /* TTX no telem -message only for first 30 seconds */
        CODELEARNING    29      /* Learning */
        CODEARTIFACT1   31      /* Level 1 Artifact noise */
        CODEARTIFACT2   32      /* same as CODESHUT? LEVEL 2 Artifact 20/30 seconds
                                   noise */

    /* End Telemetry Specific */

        CODESHUT        32      /* Shut down - all channel noise */
        CODEON           33      /* Arrhythmia is On */
        CODEOFF          34      /* Arrhythmia is Off */
        messages[2].msg_index reserved
value (short)          minutes_of_alarms_suspend

```

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Update values are sent in the Parameter Update and Extended Parameter Update structures. Values sent are expired and inspired gases, and the respiration rate value. Limits are sent in the Setup and Limits and More Setup structures. Limit alarms are sent in the Parameter update structure and messages are in the Messages structure.

CO₂ gas values and limits are stored and transmitted to the network expressed in units of mmHg. The O₂ gas value and limits are stored and transmitted expressed as a percentage, in units of 0.1. For example, a transmitted gas value of 100 would be interpreted and displayed as 10.0%. Respiration rate is expressed as breaths/minute.

If a unit conversion is necessary the barometric pressure can be found in flag [1] of the Setup and Limits structure. The barometric pressure is expressed as a byte value that is added to 530 to produce the barometric pressure in mmHg. The unit type sent in flag [0] (mmHg, %, kPa) only represents the units of the present displayed value and does not represent the units of the stored/transmitted values.

Gas values are obtained from one of a number of devices. The device is defined in flag [1] of the Setup and Limits structure. Note that O₂ is not available on all devices. For example, the sidestream and mainstream CO₂ module with Pyron or Novamatrix sensor doesn't have an O₂ parameter. Engstrom doesn't have an expired O₂ parameter or respiration rate. The codes MISSING (-32767, 0x8001) or INVALID (-32768, 0x8000) are sent if the parameter is inactive.

The gas device type can be found in the Extended Parameter Update structure. The old ID codes are still supported in the Setup and Limits structure, however if the CO₂_EXTENDED_ID bits are set the ID or device type can be found in Extended Parameter Update.

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Ventilator Parameters

Ventilator data is transferred within multiple parameter packets for the Dash, Solar 8000M/i, Solar 9500 and Unity ID.

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Update values and identification codes for the gas values are sent in the Parameter Update and Extended Parameter Update structures. Values are sent as sets of expired and inspired gases, and are identified by the expired subcode in each structure. An inactive gas value is identified by the INVALID_SUBCODE parameter in the structure.

Each of the three parameter structures may contain up to two expired/inspired sets for support of a maximum of six expired/inspired gas value sets. Under some conditions, some expired/inspired sets may be duplicated; this is signified by the high bit of the expired subcode being set.

All gas values and limits are stored and transmitted expressed as a percentage, in units of 0.01. For example, a transmitted gas value of 1000 would be interpreted and displayed as 10.00%.

If a unit conversion is necessary the barometric pressure is stored in flag[1] of the Setup and Limits structure. The barometric pressure is expressed as a byte value that is added to 530 to produce the barometric pressure in mmHg. The unit type sent in flag[0] (mmHg, %, kPa) only represents the units of the present displayed value and does not represent the units of the stored/transmitted values.

Gas values are obtained from one of a number of devices. The gas device type can be found in the Extended Parameter Update structure.

All existing devices are listed in the Extended Parameter Update structure (DATEX, EAS, MGIR, Sidestream). New devices will have the DATEX ID set in the Setup and Limits structure and the new device type code included in the Extended Parameter Update structure.

Two sets (inspired/expired) of gas limits are sent with each parcode. The Setup and Limits structure contains 1 and 1/2 sets of the limits and the More Setup structure contains 1/2 set. The subcodes for the two sets are also sent in the Setup and Limits structure and More Setup structure.

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Arterial Blood Gas Parameter

PCO₂, PO₂, pH, and HCO₃ values are transmitted under parameter code 125.

The Arterial Blood Gas values are valid when the year element of the time_stamp element of the Extended Parameter Update is not equal to 0. The completion of an ABG measurement will correspond with a change in the time_stamp value.

PCO₂ and PO₂ values are transmitted in units of mmHg. The pH value is transmitted as the actual pH value multiplied by 100. For example, a transmitted pH value of 740 would be interpreted and displayed as 7.40.

The PCO₂, PO₂, pH, or HCO₃ values may have a value of MISSING (-32767, 0x8001) or INVALID (-32768, 0x8000) if an ABG measurement has not yet been completed or a measurement error has occurred (i.e. poor pH signal quality as indicated in the par_status of the Parameter Update structure would result in the pH value being set to INVALID).

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Transcutaneous CO₂ / O₂ (Interfaced) Parameter

The values derived from the interfaces to various transcutaneous CO₂ and O₂ monitors are transmitted with the parcode tco2_par. The values, alarms, and messages from these interfaces varies with the capabilities of the device interfaced.

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Continuous Cardiac Output (Interfaced) Parameter

The values derived from the interfaces to various Continuous Cardiac Output monitors are transmitted with the parcode cco_par. The values, alarms, and messages from these interfaces varies with the capabilities of the device interfaced.

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Pulse Oximeter (Interfaced) Parameter

This structure contains the data grouped under Pulse Oximetry for external, interfaced devices, i.e., GE Medical Systems *Information Technologies* monitors. The parcodes are chosen on a first-come, first-selected basis and carry no other information. Site information is not available at this time. It is designed to map to the existing SpO₂ structures as closely as possible.

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ECG (Interfaced) Parameter

This structure contains the data grouped under ECG for external, interfaced devices, i.e., GE Medical Systems *Information Technologies* monitors. It is designed to map to the structure for hr_par as closely as possible.

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Blood Pressure (Interfaced) Parameter

This structure contains the data grouped under Invasive Blood Pressure for external, interfaced devices, i.e., GE Medical Systems *Information Technologies* monitors. The parcodes are chosen on a first-come, first-selected basis and carry no other information. The blood pressure site, e.g., ART or PA, may be determined by either the par_type field or the site value in flag[1] of the Setup and Limits structure.

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Temperature (Interfaced) Parameter

This structure contains the data grouped under Temperature for external, interfaced devices, i.e., non-GE Medical Systems *Information Technologies* monitors. The parcodes are chosen on a first-come, first-selected basis and carry no other information. Site information is not available at this time. It is designed to map to our existing temperature structures as closely as possible.

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NBP (Interfaced) Parameter

This structure contains the data grouped under Non-Invasive Blood Pressure for external, interfaced devices, i.e., GE Medical Systems *Information Technologies* monitors. It is designed to map to the existing NBP structures as closely as possible.

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Respiration (Interfaced) Parameter

This structure contains the data grouped under (impedance-based) Respiration for external, interfaced devices, i.e., GE Medical Systems *Information Technologies* monitors. It is designed to map to the `rr_par` structures as closely as possible.

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Blood Temperature/Cardiac Output (Interfaced) Parameter

This structure contains the data grouped under (bolus) Blood Temperature/Cardiac Output for external, interfaced devices, i.e., GE Medical Systems *Information Technologies* monitors. It is designed to map to the bt_par structures as closely as possible.

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Respiratory Mechanics Parameters

The Respiratory Mechanics Module provides measured ventilation data under two parcodes; rm_par and rm_par1. Some of the values are broken out to show the spontaneous (suffix _s), the mechanical (suffix _m), and the total (no suffix) components.

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More Setup

```
par_func_code (char)  PAR_MORE_SETUP_FC (3)
parcode (char)        rm_par (90)
value (short [4])
    value[0] (16 bits)
        0-2    reserved
        3-7    model, e.g., device_nova_flow_oem
        8-C    device type, device_resp_mech
        D-F    reserved
    value[1] (16 bits)
        0-7    inspired agent compensation setting (0.1%)
        8-F    expired agent compensation setting (0.1%)
    value[2] (16 bits)
        0-7    reserved
        8      rm pressure waveform available
        9      rm flow waveform available
        A      rm volume waveform available
        B-F    reserved
    value[3] (16 bits)
        0-3    sensor type (none=0, neo=1, adult=2)
        4-7    balance gas (Air=0, N2O/O2=1, O2>60%=2, HELIOX=3)
        8-F    mechanical threshold (0.5 cm H2O)
```

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SvO₂ (Interfaced) Parameter

This structure contains the data grouped under SvO₂ for external, interfaced devices, i.e., GE Medical Systems *Information Technologies* monitors. It is designed to map to the svo2_par structures as closely as possible.

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ICG Parameter

The ICG data is available when using the ICG module with DASH patient monitors. This data is transferred within three parameter groups, identified with the parcodes, icg_par, icg1_par, and icg2_par.

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BIS Module Parameter

The following structure describes the data from the BIS/EEG module when used with the BIS cable. The structure for data from the BIS interface (for the Aspect A-2000) is described elsewhere. It is identified with the parcode bism_par.

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EEG Parameter

The following structure describes the data from the BIS/EEG module when used with the EEG cable. This data is transferred within six parameter groups, identified with the parcodes, eeg_par, eeg1_par to ecg5_par.

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Information Technologies, Inc.
Information Technologies GmbH
Information Technologies Asia; GE (China) Co., Ltd.